



**Medtronic**

# IPC<sup>®</sup> System

INTEGRATED POWER CONSOLE



# Innovation *that Transforms*

Inspired by Surgeons and Nurses in the Field

The Integrated Power Console is the **ONLY** powered ENT system to offer:

- The widest range of application-specific products
- IntelliFlow® remote control for irrigation
- 0.5 ml/minute irrigation flow rates for working in the facial recess
- Integrated facial nerve monitoring with the Visao® Otologic Drill and Stim Bur Guard\*
- High-speed, application-specific curved burs for otologic procedures
- A complete anterior skull base surgery solution

\*For use with Medtronic's NIM® 2.0 and 3.0 nerve monitoring systems

# Integrated Power Console (IPC®) System *Innovati*

The new Integrated Power Console delivers innovation that transforms your O.R. Its design was inspired by feedback from you, the surgeons and nurses in the field. Our IPC® system is the only powered ENT surgery console with the most expansive range of application-specific products. From inferior turbinoplasties to complex skull base cases, the IPC platform is your ENT surgery system.

*Integrated facial nerve monitoring with the Visao® Otologic Drill and Stim Bur Guard\**

*Easy-to-use touchscreen menu for each handpiece*

*IntelliFlow® remote control for irrigation*

*Maximum visualization for anterior skull base*

*Connect multiple handpieces at once*



\*For use with Medtronic's NIM® 2.0 and 3.0 nerve monitoring systems

ion that Transforms

## Complete Integration

*Power all of your favorite Medtronic ENT handpieces*

- Straightshot® M4 Microdebrider
- Endo-Scrub® 2 Lens Cleaning Sheaths
- Skeeter® Otologic Drill
- Visao® High-Speed Otologic Drill



## Flexible Irrigation

*On-drill irrigation reinvented*

- Disposable IntelliFlow® remote allows precise flow control from sterile field
- 0.5 ml/minute irrigation flow rates for working in the facial recess
- Snap-on cassette for easier set-up
- Inner diameter of tubing enables continuous flow, even at low flow rates



## Multifunction Footpedal

*Ergonomic control at your feet*

- Three buttons provide maximum control
- Built-in footrests and ergonomic angle for comfort
- Hoop locks firmly in place, then collapses easily for storage
- Non-slip design and compact size is ideal for your O.R.



# Indigo™

## High Speed Otologic Drill

- High speed and ample torque for fast debulking of the mastoid
- Available with straight and angled attachments
- Improved surgical access and visibility in narrow otologic spaces
- Works with a wide variety of Medtronic Xomed otology burs
- Grooves on attachment nose design allows for more stable irrigation



## Anterior Skull Base Surgery

### *The complete solution for superior visualization*


- Curved burs and blades allow better visibility and access
- Integrated suction keeps surgical site clear

*The IPC platform is the ONLY powered ENT system that provides the widest range of application-specific products.*

## Stim Bur Guard for the Visao® Otologic Drill

### *Integrated facial nerve monitoring*

- Provides 1-3 mm advance warning of facial nerve proximity
- Helps reduce the risk of patient injury for added peace-of-mind
- Works with our NIM® intraoperative nerve monitor\*



\*For use with Medtronic's NIM® 2.0 and 3.0 nerve monitoring systems

## Ordering Information

### IPC® System and Endo-Scrub® 2 Accessories

Product	Description	Qty
1898001	Integrated Power Console	1
1898430	Integrated Power Console Multifunction Footpedal	1
1897510	Basket	1
1852000	Footswitch Endo-Scrub 2	1
1991015	Finger Switch Endo-Scrub 2	1

### IPC® System Power Cords

Product	Description	Qty
1895820	Power Cord, Standard, North America, 3M	1
1895821	Power Cord, Standard, North America, 6M	1
1895824	Power Cord, UK/Ireland, 6M	1
1895825	Power Cord, Conintental Europe, 6M	1
1895823	Power Cord, Japan	1

### IPC® System Manuals

Product	Description	Qty
1898851	IPC Manual (English Only)	1
1898851A	IPC Manual (EL, EN, ES, FR, PT)	1
1898851B	IPC Manual (DE, EN, FR, IT, NL)	1
1898851C	IPC Manual (DA, EN, FI, NO, SV)	1
1898851D	IPC Manual (CS, EN, HU, PL, TR)	1
1898851E	IPC Manual (EN, RU) - CD only	1

### Straightshot® M4 Microdebrider Components

Product	Description	Qty
1898200T	XPS® Straightshot® M4 High-Speed Microdebrider	1
1898400	XPS Straightshot M4 Instrument Tray	1

For further information, please call Medtronic ENT at 800-874-5797 or 904-296-9600. You may also consult our website at [www.MedtronicENT.com](http://www.MedtronicENT.com).

#### Medtronic ENT

Medtronic USA, Inc.  
6743 Southpoint Drive North  
Jacksonville, FL 32216 USA  
[www.MedtronicENT.com](http://www.MedtronicENT.com)  
Telephone: (904) 296-9600  
Toll free: (800) 874-5797  
Fax: (904) 296-9666

#### International Telephone Numbers

Australia 1-800-668-670	Korea 82-2-3404-3600
Belgium 32-2456-09-00	Lebanon 961-1-370-670
Canada 1-800-217-1617	Luxembourg 32-2456-09-09
China 86-21-50800998	Netherlands 31-45-566-8000
France 33-470-679-800	Poland 48-22-465-6942
Germany 49-211-5293-209	Singapore 65-6776-6255
Hong Kong 852-2919-1312	Spain 34-91-625-05-40
India 91-22-26836733	UK 44-1923-212-213
Italy 39-02-24137-324	USA 1-904-296-9600
Japan 81-6-4795-1506	

### Visao® High-Speed Otologic Drill Components

Product	Description	Qty
3334800	Visao High-Speed Otologic Drill	1
3334625	Visao High-Speed Otologic Drill Extended Bur Guard without Irrigation	1
3334635	Visao High-Speed Otologic Drill Extended Bur Guard with Irrigation	1
3334895	Visao High-Speed Otologic Drill Instrument Tray	1
3155617	Bur Brush	1
3112500	Micro-Craft Chuck Brush Cleaner	60
3318603	Visao Irrigation Tubing with IntelliFlow Remote Control	1
3318605	Suction/Irrigation Connector Set	5
3318606	Visao Coolant Tubing	2
3318601	Stim Bur Guard with Irrigation	1
3318602	Stim Bur Guard without Irrigation	1

### Skeeter® Components

Product	Description	Qty
3055601	Skeeter Ultralite Handpiece	1
3155615	Skeeter Instrument Tray	1
3052460	Skeeter Otologic Metric Scale	1
3155641	Skeeter Oto-Flex Bur, Orange/White, Diamond, 1 mm	1
3155631	Skeeter Oto-Flex Bur, Orange, Carbide, 1 mm	1
3155647	Skeeter Oto-Flex Bur, Green/White, Diamond, 0.7 mm	1
3155632	Skeeter Oto-Flex Bur, Red, Carbide, 2.3 mm	1
3155648	Skeeter Oto-Flex Bur, Yellow/White, Diamond, 0.8 mm	1
3155642	Skeeter Oto-Flex Bur, Red/White, Diamond, 2.3 mm	1
3155646	Skeeter Oto-Flex Bur, Blue/White, Diamond, 0.6 mm	1
3155643	Skeeter Oto-Flex Bur, Gray/White, Diamond, 1.4 mm	1
3155645	Skeeter Oto-Flex Bur, Brown/White, Diamond, 1.8 mm	1
3155635	Skeeter Oto-Flex Bur, Violet, Carbide, 0.5 mm	1

### Indigo™ Components

Product	Description	Qty
1845000	Indigo High-Speed Otologic Drill	1
1845010	Indigo Otologic Straight Attachment	1
1845020	Indigo Otologic Angled Attachment	1
1845040	Indigo Instrument Tray	1
1845030	Indigo Tubing	1

This literature is intended for the exclusive use of physicians. Rx only.  
©2013, Medtronic, Inc. \*Registered mark of Medtronic, Inc.  
UC201400184 EN 06.13

# ENT MR8™ HIGH SPEED DRILL CATALOG



# ENT MR8™ High Speed Drill System

## *Table of contents*

<b>Introduction</b> .....	<b>4</b>
<b>Bur Offering</b> .....	<b>5</b>
Cutting Fine.....	6
Cutting Standard.....	7
Cutting Carbide .....	8
Diamond Fine .....	9
Diamond Medium.....	10
Diamond Coarse .....	11
Diamond Extra Coarse.....	11
ENT MR8™ Midas Rex ClearView™ Curved Burs .....	12
Curved Burs Cutting .....	13
Curved Burs Fine Diamond.....	14
Curved Burs Coarse Diamond.....	14
<b>Accessories</b> .....	<b>15</b>

# WORLD-CLASS TECHNOLOGY IN YOUR HANDS



The right high speed drill helps you do your best work during otolaryngology procedures.

The new Medtronic ENT MR8™ drill features the latest Midas Rex™ technology to enable precise control and improve visualization. A generational advancement,\* it offers the best of previous Medtronic drills\* — with a focus on comfort, efficiency and responsiveness — all on a modern platform.

## ORDERING INFORMATION

ENT MR8 CFN	ENT MR8 DESCRIPTION
1847DRLMTR	DRILL MOTOR ENTMR8 1847DRLMTR OTOLGY
1847ATTSS	ATTACHMENT ENTMR8 1847ATTSS STRT STD
1847ATTAS	ATTACHMENT ENTMR8 1847ATTAS ANG STD
1847ATTSV	ATTACHMENT ENTMR8 1847ATTSV STRT VAR
1847ATTAV	ATTACHMENT ENTMR8 1847ATTAV ANG VAR
1847IRTUBBX	TUBING ENTMR8 1847IRTUBEBX IRRIGATION 5P
1847IRTUBEA	TUBING ENTMR8 1847IRTUBEEA IRRIGATION EA
1847CVTUBBX	TUBING ENTMR8 1847CVTUBBX CLEARVIEW 5PK
1847CVTUBEA	TUBING ENTMR8 1847CVTUBEA CLEARVIEW EA

\* When compared to Legend™, Indigo™, and Visao™ High Speed Drills.

# ENT MR8™ High Speed Drill System

*The Power Behind ENT MR8™*

1

## Motor

- Single-click Attachment lock
- Better handpiece balance
- Continuously drilling with a rotational speed of up to 75,000 RPM



4

## Attachments

- Fixed straight and angled configurations
- Expanded offering include variable exposure attachments.
  - These attachments allow you to vary dissecting tool exposure.
  - They can increase exposure up to 8mm.
- Streamlined attachment profiles to maximize visibility



94

## Dissecting Burs

- Expanded offering to include fine cutters and extra coarse diamond
- Single Use
- Disposable bur line not backward-compatible with Indigo™/Visao™ system



28

## ENT MR8™ Midas Rex Clearview™ Curved Burs

- Self-irrigating attachment-tool solution
- Extended self-irrigating sleeve
- Direct connect burs
- Single use





## New Accessories


- Optional sterilization container with filter

# ENT MR8™ High Speed Drill

## Types of Burs



## BUR CFN NOMENCLATURE



### Family

B	Regular Bur
CV	ENT MR8™ Midas Rex ClearView™ Bur

### Bur Tip Size

0.5 – 8.0 mm	Diameter in Millimeters (60=6.0 mm, 25=2.5 mm)
--------------	---

### Bur Length Code

S	Standard
L	Long

### Bur Type and Coarse Description

CT	Cutting (Standard)
CF	Cutting Fine
CC	Cutting Carbide
DF	Diamond Fine
DM	Diamond Medium
DC	Diamond Coarse
DX	Diamond Extra Coarse

# ENT MR8™

## Cutting Burs



### Fine

Fine Cutting Bur's less aggressive grind features optimized flute geometry. Fine Cutting Burs are designed for slower bone removal with increased control.

ENT MR8 CFN	ENT MR8 DESCRIPTION	TIP DIAMETER (mm)	TIP MATERIAL	LENGTH
1847B60SCF	BUR ENTMR8 1847B60SCF CUT FINE 6MM STD	6	STAINLESS STEEL	STANDARD
1847B50SCF	BUR ENTMR8 1847B50SCF CUT FINE 5MM STD	5	STAINLESS STEEL	STANDARD
1847B40SCF	BUR ENTMR8 1847B40SCF CUT FINE 4MM STD	4	STAINLESS STEEL	STANDARD
1847B30SCF	BUR ENTMR8 1847B30SCF CUT FINE 3MM STD	3	STAINLESS STEEL	STANDARD
1847B25SCF	BUR ENTMR8 1847B25SCF CUT FINE 2.5MM STD	2.5	STAINLESS STEEL	STANDARD
1847B20SCF	BUR ENTMR8 1847B20SCF CUT FINE 2MM STD	2	STAINLESS STEEL	STANDARD
1847B15SCF	BUR ENTMR8 1847B15SCF CUT FINE 1.5MM STD	1.5	STAINLESS STEEL	STANDARD
1847B10SCF	BUR ENTMR8 1847B10SCF CUT FINE 1MM STD	1	STAINLESS STEEL	STANDARD
1847B35LCF	BUR ENTMR8 1847B35LCF CUT FINE 3.5MM LNG	3.5	TOOL STEEL	LONG
1847B30LCF	BUR ENTMR8 1847B30LCF CUT FINE 3MM LNG	3	TOOL STEEL	LONG
1847B25LCF	BUR ENTMR8 1847B25LCF CUT FINE 2.5MM LNG	2.5	TOOL STEEL	LONG
1847B20LCF	BUR ENTMR8 1847B20LCF CUT FINE 2MM LNG	2	TOOL STEEL	LONG
1847B15LCF	BUR ENTMR8 1847B15LCF CUT FINE 1.5MM LNG	1.5	STAINLESS STEEL	LONG
1847B10LCF	BUR ENTMR8 1847B10LCF CUT FINE 1MM LNG	1	STAINLESS STEEL	LONG

# ENT MR8™

## Cutting Burs



### Standard

Standard Cutting bur's speed and aggressiveness will be comparable to previous Medtronic cutting burs. Tool steel's increased hardness will allow better edge and sharpness retention throughout procedures when compared to stainless steel burs.

ENT MR8 CFN	ENT MR8 DESCRIPTION	TIP DIAMETER (mm)	TIP MATERIAL	LENGTH
1847B80SCT	BUR ENTMR8 1847B80SCT CUT 8MM STD	8	STAINLESS STEEL	STANDARD
1847B70SCT	BUR ENTMR8 1847B70SCT CUT 7MM STD	7	STAINLESS STEEL	STANDARD
1847B60SCT	BUR ENTMR8 1847B60SCT CUT 6MM STD	6	TOOL STEEL	STANDARD
1847B50SCT	BUR ENTMR8 1847B50SCT CUT 5MM STD	5	TOOL STEEL	STANDARD
1847B40SCT	BUR ENTMR8 1847B40SCT CUT 4MM STD	4	TOOL STEEL	STANDARD
1847B35SCT	BUR ENTMR8 1847B35SCT CUT 3.5MM STD	3.5	TOOL STEEL	STANDARD
1847B30SCT	BUR ENTMR8 1847B30SCT CUT 3MM STD	3	TOOL STEEL	STANDARD
1847B25SCT	BUR ENTMR8 1847B25SCT CUT 2.5MM STD	2.5	TOOL STEEL	STANDARD
1847B20SCT	BUR ENTMR8 1847B20SCT CUT 2MM STD	2	TOOL STEEL	STANDARD
1847B15SCT	BUR ENTMR8 1847B15SCT CUT 1.5MM STD	1.5	TOOL STEEL	STANDARD
1847B10SCT	BUR ENTMR8 1847B10SCT CUT 1MM STD	1	TOOL STEEL	STANDARD
1847B35LCT	BUR ENTMR8 1847B35LCT CUT 3.5MM LNG	3.5	TOOL STEEL	LONG
1847B30LCT	BUR ENTMR8 1847B30LCT CUT 3MM LNG	3	TOOL STEEL	LONG
1847B25LCT	BUR ENTMR8 1847B25LCT CUT 2.5MM LNG	2.5	TOOL STEEL	LONG
1847B20LCT	BUR ENTMR8 1847B20LCT CUT 2MM LNG	2	TOOL STEEL	LONG
1847B15LCT	BUR ENTMR8 1847B15LCT CUT 1.5MM LNG	1.5	TOOL STEEL	LONG
1847B10LCT	BUR ENTMR8 1847B10LCT CUT 1MM LNG	1	TOOL STEEL	LONG
1847B05LCT	BUR ENTMR8 1847B05LCT CUT 0.5MM LNG	0.5	TOOL STEEL	LONG

# ENT MR8™

## Cutting Burs



### Carbide

Carbide Cutting Burs have a similar grind to the standard cutting burs. Carbide Cutting Burs will have the best edge and sharpness retention for the most detailed or demanding procedures.

ENT MR8 CFN	ENT MR8 DESCRIPTION	TIP DIAMETER (mm)	TIP MATERIAL	LENGTH
1847B80SCC	BUR ENTMR8 1847B80SCC CUT CBD 8MM STD	8	CARBIDE	STANDARD
1847B70SCC	BUR ENTMR8 1847B70SCC CUT CBD 7MM STD	7	CARBIDE	STANDARD
1847B60SCC	BUR ENTMR8 1847B60SCC CUT CBD 6MM STD	6	CARBIDE	STANDARD
1847B50SCC	BUR ENTMR8 1847B50SCC CUT CBD 5MM STD	5	CARBIDE	STANDARD
1847B40SCC	BUR ENTMR8 1847B40SCC CUT CBD 4MM STD	4	CARBIDE	STANDARD
1847B30SCC	BUR ENTMR8 1847B30SCC CUT CBD 3MM STD	3	CARBIDE	STANDARD
1847B20SCC	BUR ENTMR8 1847B20SCC CUT CBD 2MM STD	2	CARBIDE	STANDARD
1847B10SCC	BUR ENTMR8 1847B10SCC CUT CBD 1MM STD	1	CARBIDE	STANDARD

# ENT MR8™

## Diamond Burs



Fine

ENT MR8 CFN	ENT MR8 DESCRIPTION	TIP DIAMETER (mm)	LENGTH
1847B60SDF	BUR ENTMR8 1847B60SDF DMD FINE 6MM STD	6	STANDARD
1847B50SDF	BUR ENTMR8 1847B50SDF DMD FINE 5MM STD	5	STANDARD
1847B40SDF	BUR ENTMR8 1847B40SDF DMD FINE 4MM STD	4	STANDARD
1847B30SDF	BUR ENTMR8 1847B30SDF DMD FINE 3MM STD	3	STANDARD
1847B25SDF	BUR ENTMR8 1847B25SDF DMD FINE 2.5MM STD	2.5	STANDARD
1847B20SDF	BUR ENTMR8 1847B20SDF DMD FINE 2MM STD	2	STANDARD
1847B15SDF	BUR ENTMR8 1847B15SDF DMD FINE 1.5MM STD	1.5	STANDARD
1847B10SDF	BUR ENTMR8 1847B10SDF DMD FINE 1MM STD	1	STANDARD
1847B50LDF	BUR ENTMR8 1847B50LDF DMD FINE 5MM LNG	5	LONG
1847B40LDF	BUR ENTMR8 1847B40LDF DMD FINE 4MM LNG	4	LONG
1847B30LDF	BUR ENTMR8 1847B30LDF DMD FINE 3MM LNG	3	LONG
1847B25LDF	BUR ENTMR8 1847B25LDF DMD FINE 2.5MM LNG	2.5	LONG
1847B20LDF	BUR ENTMR8 1847B20LDF DMD FINE 2MM LNG	2	LONG
1847B15LDF	BUR ENTMR8 1847B15LDF DMD FINE 1.5MM LNG	1.5	LONG
1847B10LDF	BUR ENTMR8 1847B10LDF DMD FINE 1MM LNG	1	LONG

# ENT MR8™

## Diamond Burs



Medium

ENT MR8 CFN	ENT MR8 DESCRIPTION	TIP DIAMETER (mm)	LENGTH
1847B80SDM	BUR ENTMR8 1847B80SDM DMD MED 8MM STD	8	STANDARD
1847B70SDM	BUR ENTMR8 1847B70SDM DMD MED 7MM STD	7	STANDARD
1847B60SDM	BUR ENTMR8 1847B60SDM DMD MED 6MM STD	6	STANDARD
1847B50SDM	BUR ENTMR8 1847B50SDM DMD MED 5MM STD	5	STANDARD
1847B40SDM	BUR ENTMR8 1847B40SDM DMD MED 4MM STD	4	STANDARD
1847B35SDM	BUR ENTMR8 1847B35SDM DMD MED 3.5MM STD	3.5	STANDARD
1847B30SDM	BUR ENTMR8 1847B30SDM DMD MED 3MM STD	3	STANDARD
1847B25SDM	BUR ENTMR8 1847B25SDM DMD MED 2.5MM STD	2.5	STANDARD
1847B20SDM	BUR ENTMR8 1847B20SDM DMD MED 2MM STD	2	STANDARD
1847B15SDM	BUR ENTMR8 1847B15SDM DMD MED 1.5MM STD	1.5	STANDARD
1847B10SDM	BUR ENTMR8 1847B10SDM DMD MED 1MM STD	1	STANDARD
1847B40LDM	BUR ENTMR8 1847B40LDM DMD MED 4MM LNG	4	LONG
1847B35LDM	BUR ENTMR8 1847B35LDM DMD MED 3.5MM LNG	3.5	LONG
1847B30LDM	BUR ENTMR8 1847B30LDM DMD MED 3MM LNG	3	LONG
1847B25LDM	BUR ENTMR8 1847B25LDM DMD MED 2.5MM LNG	2.5	LONG
1847B20LDM	BUR ENTMR8 1847B20LDM DMD MED 2MM LNG	2	LONG
1847B15LDM	BUR ENTMR8 1847B15LDM DMD MED 1.5MM LNG	1.5	LONG
1847B12LDM	BUR ENTMR8 1847B12LDM DMD MED 1.2MM LNG	1.2	LONG
1847B10LDM	BUR ENTMR8 1847B10LDM DMD MED 1MM LNG	1	LONG
1847B08LDM	BUR ENTMR8 1847B08LDM DMD MED 0.8MM LNG	0.8	LONG
1847B06LDM	BUR ENTMR8 1847B06LDM DMD MED 0.6MM LNG	0.6	LONG

# ENT MR8™

## Diamond Burs



### Coarse

ENT MR8 CFN	ENT MR8 DESCRIPTION	TIP DIAMETER (mm)	LENGTH
1847B70SDC	BUR ENTMR8 1847B70SDC DMD CRS 7MM STD	7	STANDARD
1847B60SDC	BUR ENTMR8 1847B60SDC DMD CRS 6MM STD	6	STANDARD
1847B50SDC	BUR ENTMR8 1847B50SDC DMD CRS 5MM STD	5	STANDARD
1847B40SDC	BUR ENTMR8 1847B40SDC DMD CRS 4MM STD	4	STANDARD
1847B30SDC	BUR ENTMR8 1847B30SDC DMD CRS 3MM STD	3	STANDARD
1847B25SDC	BUR ENTMR8 1847B25SDC DMD CRS 2.5MM STD	2.5	STANDARD
1847B20SDC	BUR ENTMR8 1847B20SDC DMD CRS 2MM STD	2	STANDARD
1847B15SDC	BUR ENTMR8 1847B15SDC DMD CRS 1.5MM STD	1.5	STANDARD
1847B10SDC	BUR ENTMR8 1847B10SDC DMD CRS 1MM STD	1	STANDARD
1847B40LDC	BUR ENTMR8 1847B40LDC DMD CRS 4MM LNG	4	LONG
1847B35LDC	BUR ENTMR8 1847B35LDC DMD CRS 3.5MM LNG	3.5	LONG
1847B30LDC	BUR ENTMR8 1847B30LDC DMD CRS 3MM LNG	3	LONG
1847B20LDC	BUR ENTMR8 1847B20LDC DMD CRS 2MM LNG	2	LONG
1847B15LDC	BUR ENTMR8 1847B15LDC DMD CRS 1.5MM LNG	1.5	LONG
1847B10LDC	BUR ENTMR8 1847B10LDC DMD CRS 1MM LNG	1	LONG




### Extra Coarse

ENT MR8 CFN	ENT MR8 DESCRIPTION	TIP DIAMETER (mm)	LENGTH
1847B60SDX	BUR ENTMR8 1847B60SDX DMD XCRS 6MM STD	6	STANDARD
1847B50SDX	BUR ENTMR8 1847B50SDX DMD XCRS 5MM STD	5	STANDARD
1847B40SDX	BUR ENTMR8 1847B40SDX DMD XCRS 4MM STD	4	STANDARD

# ENT MR8™ Midas Rex ClearView™ Curved Burs

ENT MR8™ Midas Rex ClearView™ curved burs feature responsive, integrated self-irrigating channels close to the cutting tip to irrigate the surgical site without impeding the view.



# ENT MR8™ Midas Rex ClearView™ Curved Burs



## Cutting

ENT MR8 CFN	ENT MR8 DESCRIPTION	TIP DIAMETER (mm)	TIP MATERIAL	LENGTH
1847CV10SCT	BUR ENTMR8 1847CV10SCT CUT 1MM STD	1	TOOL STEEL	STANDARD
1847CV15SCT	BUR ENTMR8 1847CV15SCT CUT 1.5MM STD	1.5	TOOL STEEL	STANDARD
1847CV20SCT	BUR ENTMR8 1847CV20SCT CUT 2MM STD	2	TOOL STEEL	STANDARD
1847CV30SCT	BUR ENTMR8 1847CV30SCT CUT 3MM STD	3	TOOL STEEL	STANDARD
1847CV40SCT	BUR ENTMR8 1847CV40SCT CUT 4MM STD	4	TOOL STEEL	STANDARD
1847CV10LCT	BUR ENTMR8 1847CV10LCT CUT 1MM LNG	1	TOOL STEEL	LONG
1847CV15LCT	BUR ENTMR8 1847CV15LCT CUT 1.5MM LNG	1.5	TOOL STEEL	LONG
1847CV20LCT	BUR ENTMR8 1847CV20LCT CUT 2MM LNG	2	TOOL STEEL	LONG
1847CV30LCT	BUR ENTMR8 1847CV30LCT CUT 3MM LNG	3	TOOL STEEL	LONG
1847CV40LCT	BUR ENTMR8 1847CV40LCT CUT 4MM LNG	4	TOOL STEEL	LONG

Irrigation tubing sets to be used with the ENT MR8™ Midas Rex Clearview™ Burs:

ENT MR8 CFN	ENT MR8 DESCRIPTION
1847CVTUBBX	TUBING ENTMR8 1847CVTUBBX CLEARVIEW 5PK
1847CVTUBEA	TUBING ENTMR8 1847CVTUBEA CLEARVIEW EA

# ENT MR8™ Midas Rex ClearView™

## Curved Burs



### Fine Diamond

ENT MR8 CFN	ENT MR8 DESCRIPTION	TIP DIAMETER (mm)	LENGTH
1847CV05SDF	BUR ENTMR8 1847CV05SDF DMD 0.5MM STD	0.5	STANDARD
1847CV08SDF	BUR ENTMR8 1847CV08SDF DMD 0.8MM STD	0.8	STANDARD
1847CV10SDF	BUR ENTMR8 1847CV10SDF DMD 1MM STD	1	STANDARD
1847CV12SDF	BUR ENTMR8 1847CV12SDF DMD 1.2MM STD	1.2	STANDARD
1847CV15SDF	BUR ENTMR8 1847CV15SDF DMD 1.5MM STD	1.5	STANDARD
1847CV20SDF	BUR ENTMR8 1847CV20SDF DMD 2MM STD	2	STANDARD
1847CV30SDF	BUR ENTMR8 1847CV30SDF DMD 3MM STD	3	STANDARD
1847CV40SDF	BUR ENTMR8 1847CV40SDF DMD 4MM STD	4	STANDARD
1847CV10LDF	BUR ENTMR8 1847CV10LDF DMD 1MM LNG	1	LONG
1847CV20LDF	BUR ENTMR8 1847CV20LDF DMD 2MM LNG	2	LONG
1847CV30LDF	BUR ENTMR8 1847CV30LDF DMD 3MM LNG	3	LONG
1847CV40LDF	BUR ENTMR8 1847CV40LDF DMD 4MM LNG	4	LONG



### Coarse Diamond

ENT MR8 CFN	ENT MR8 DESCRIPTION	TIP DIAMETER (mm)	LENGTH
1847CV20SDC	BUR ENTMR8 1847CV20SDC DMD 2MM STD	2	STANDARD
1847CV30SDC	BUR ENTMR8 1847CV30SDC DMD 3MM STD	3	STANDARD
1847CV40SDC	BUR ENTMR8 1847CV40SDC DMD 4MM STD	4	STANDARD
1847CV20LDC	BUR ENTMR8 1847CV20LDC DMD 2MM LNG	2	LONG
1847CV30LDC	BUR ENTMR8 1847CV30LDC DMD 3MM LNG	3	LONG
1847CV40LDC	BUR ENTMR8 1847CV40LDC DMD 4MM LNG	4	LONG

# Accessories

CFN	Description
1847BRUSH	BRUSH ENTMR8 1847BRUSH CLEANING 2.4MM 5P
1847TRAY	TRAY ENTMR8 1847TRAY INSTRUMENT 1/2 DIN
1847CNTR	CONTAINER ENTMR8 1847CNTR STERIL 1/2 DIN
1847BRACK	RACK ENTMR8 1847BRACK BUR RACK
1847WRENCH	WRENCH ENTMR8 1847WRENCH MOTOR
1847IRTUBBX	TUBING ENTMR8 1847IRTUBBX IRRIGATION 5PK
1847IRTUBEA	TUBING ENTMR8 1847IRTUBEA IRRIGATION EA



1847BRACK



1847IRTUBEA



1847CNTR



1847TRAY

Rx only. Refer to product instruction manual/package insert for instructions, warnings, precautions and contraindications.

For further information, please call Medtronic ENT at 800.874.5797 or consult Medtronic's website at [www.medtronicent.com](http://www.medtronicent.com).

## Medtronic

### ENT

6743 Southpoint Drive N  
Jacksonville, FL 32216  
USA  
Toll free: (800) 874-5797  
Telephone: (904) 296-9600  
Fax: (800) 678-3995

### International Telephone Numbers

Adriatic Region 385-1-488-1120	Israel 972-9-972-4400
Argentina 54-11-4898-5700	Italy 39-02-24137-324
Australia 1-800-668-670	Japan 81-3-6430-2017
Baltic Region 37-1-67560226	Korea 82-2-3404-3600
Belgium 32-2456-09-09	Lebanon 961-1-370-670
Brazil 55-11-2182-9200	Luxembourg 32-2456-09-09
Canada 1-800-217-1617	Malaysia 60-37-953-4800
Chile 56-2-2655-5110	Mexico 52-55-11-02-90-30
China 86-21-20325888	Netherlands 31-45-566-8800
Colombia 57-1-742-7300	Poland 48-22-4656900
Czech Republic 420-2-9657-9580	Russian Fed. 7-495-580-73-77
France 33-470-679-800	Singapore 65-6776-6255
Germany 49-2159-8149-209	South Africa 27-11-466-1820
Greece 30-210-67-79-099	Spain 34-91-625-05-40
Hong Kong 852-2919-1300	Taiwan 886-2-2183-6000
Hungary 36-188-90600	Thailand 662-232-7400
India 91-22-26836733	UK 44-1923-205-166

© 2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. UC202105691 EN